

# Mag-Bind® PX Blood RNA 96 Kit

Fast & reliable RNA extraction from blood samples stabilized in PAXgene™ or Tempus™ tubes

CAT NO: M7763



## 96 samples

96 samples in 1 hour after initial spin



## Quality

High quality RNA with minimal gDNA contamination



## Bead-Based

Manual or automated processing



## Automatable


Adaptable on most open-ended liquid handlers

## Cost-Effective

**30%** less than the competition on average

**Mag-Bind® PX Blood RNA 96 Kit** provides a convenient and fast RNA extraction solution for samples stabilized in PAXgene™ blood RNA tubes or Applied Biosystems' Tempus™ blood RNA tubes. This system combines Mag-Bind® particles and HiBind® silica technologies for purification of total RNA from up to 2.5 mL preserved blood samples.

The blood sample is spun down and the crude RNA/DNA pellet is collected and washed. The pellet is then resuspended and digested with Proteinase K. The lysate is applied to a filter plate to remove genomic DNA. This procedure completely removes inhibitors and sample stabilization reagents to allow for reliable downstream analysis. High quality purified RNA can be used for downstream applications such as qRT-PCR, RT-PCR, and microarray analysis.



## Electropherogram Representing Purified RNA

**Figure 1.** 2.5 mL whole blood was drawn to PAXgene™ blood RNA tubes and stored at room temperature for 24 hours. RNA was purified following manufacturer's recommended protocols. RNA was analyzed on Agilent's TapeStation® 2200.

## RNA Integrity Assessment

**Table 1.** The quality of the RNA purified was assessed by the RNA integrity score (RIN<sup>e</sup>) determined from the electrophoretic trace and Agilent's TapeStation® 2200 system.

Manufacturer	RIN <sup>e</sup>	28S/18S Ratio	Conc. (ng/μL)
Omega Bio-tek	7.5	2.1	65.3
Omega Bio-tek	7.6	2.2	29.5
Company Q	7.2	2.2	33.3
Company Q	7.1	2.2	34.7

Product Description	Preps	Cat No.
Mag-Bind® PX Blood RNA 96 Kit	1 x 96	M7763-00
	4 x 96	M7763-01

For free samples of any of our kits, visit [www.omegabiotek.com](http://www.omegabiotek.com)



innovations in nucleic acid isolation

Omega Bio-tek, Inc.  
400 Pinnacle Way, Suite 450  
Norcross, GA 30071

Phone: 770-931-8400  
Email: [info@omegabiotek.com](mailto:info@omegabiotek.com)  
Web: [www.omegabiotek.com](http://www.omegabiotek.com)



© 2018 Omega Bio-tek, Inc. All rights reserved. E.Z.N.A., E-Z 96 and Mag-Bind are trademarks of Omega Bio-tek, Inc. All other trademarks and trade names are the property of their respective holders. For research use only. Lit No. SL-0045