



E.Z.N.A.[®] Gel Extraction Kit

D2500-00	5 preps
D2500-01	50 preps
D2500-02	200 preps

Manual Date: April 2024
Revision Number: v4.1

For Research Use Only

E.Z.N.A.® Gel Extraction Kit

Table of Contents

Introduction.....	2
Illustrated Protocol.....	3
Kit Contents/Storage and Stability.....	4
Warnings and Safety Information.....	5
Preparing Reagents.....	6
Guideline for Vacuum Manifold.....	7
Spin Protocol.....	8
Vacuum Protocol.....	11
Enzymatic Reaction Protocol.....	14
Troubleshooting Guide.....	16

Manual Date: April 2024

Revision Number: v4.1



Introduction

The E.Z.N.A.® family of products is an innovative system that radically simplifies the extraction and purification of nucleic acids from a variety of sources. The key to this system is the HiBind® matrix that specifically, but reversibly, binds DNA or RNA under optimized conditions allowing proteins and other contaminants to be removed. Nucleic acids are easily eluted with deionized water or a low salt buffer.

Gel purification of DNA is a common technique for the isolation of specific fragments from reaction mixtures. However, most methods either fail to completely remove agarose or shear DNA which can lead to problems in downstream manipulations. The E.Z.N.A.® Gel Extraction Kit uses proprietary chemistry and HiBind® technology to recover DNA fragments between 70 bp and 20 kb with yields exceeding 85%. The DNA band of interest is excised from the gel, dissolved in Binding Buffer, and transferred to a HiBind® DNA Mini Column. Following three rapid wash steps, DNA is eluted with the Elution Buffer and is ready for other applications. DNA is suitable for ligations, PCR, sequencing, restriction digestion, or various labeling reactions. In addition, this kit can be also used to recover DNA directly from enzymatic reactions such as PCR and enzyme digestion reactions.

Benefits of the E.Z.N.A.® Gel Extraction Kit

- **Fast** - DNA recovery from agarose gel < 10 minutes
- **Reliability** - Optimized buffers that guarantee pure DNA
- **Safety** - No organic extractions
- **Quality** - Purified DNA is suitable for most applications

Binding Capacity

Each HiBind® DNA Mini Column can bind ~25 µg DNA.

New In this Edition

April 2024:

- Addition of Warnings and Safety Information.


January 2019:

- SPW Wash Buffer has been renamed SPW Buffer and Binding Buffer (XP2) has been renamed XP2 Binding Buffer. These are name changes only. The formulations have not changed.
- All sizes of D2501 have been discontinued.

November 2018:

- Omega Bio-tek's logo has changed.
- VAC-08 has been discontinued and is no longer available for purchase.

Spin Protocol



Excise DNA Band
from Gel



Incubate and Melt Gel
with Binding Buffer



Transfer Lysate & Bind



Wash 3x




Dry



Elute

Vacuum/Spin Protocol




Excise DNA Band
from Gel



Incubate and Melt Gel
with Binding Buffer



Transfer Lysate & Bind



Wash 3x



Dry



Elute

Kit Contents and Storage

Product	D2500-00	D2500-01	D2500-02
Purifications	5	50	200
HiBind® DNA Mini Columns	5	50	200
2 mL Collection Tubes	5	50	200
XP2 Binding Buffer	5 mL	40 mL	150 mL
SPW Buffer	5 mL	25 mL	3 x 25 mL
Elution Buffer	2 mL	30 mL	30 mL
User Manual	✓	✓	✓

Storage and Stability

All E.Z.N.A.® Gel Extraction Kit components are guaranteed for at least 12 months from the date of purchase when stored at room temperature. Please ensure that the bottle of XP2 Binding Buffer is tightly capped when not in use. Store all other components at room temperature and away from bright light. If any precipitates form in the buffers, warm at 37°C to dissolve.

Warnings and Safety Information

Warnings

This kit is for research use only.

Please read all instructions carefully before using the kit.

Decontaminate and dispose of all potentially infectious materials in accordance with applicable local, state, and national regulations. Please refer to safety data sheets (SDSs) for information on disposal of different components included in this kit.

Safety Information

All chemicals and biological materials are potentially hazardous. Biological samples such as plasma, serum, tissues, body fluids, blood etc. are potentially infectious and must be treated as biohazardous materials. Conduct all work in properly equipped facilities following universal precautions and using appropriate personal safety equipment such as disposable gloves, lab coats, safety glasses etc. as required by policies and procedures outlined by your facility. Please refer to safety data sheets (SDSs) for information on safe handling, transport and disposal of different components included in this kit. SDSs are made available in PDF format on the product page at www.omegabiotek.com. Discard all waste in accordance with the local safety regulations.

Some of the buffers included in the product contain guanidine-based chaotropic agents, which can form highly reactive compounds when combined with bleach. DO NOT add bleach or acidic solutions to guanidine-containing waste. Please access the SDSs online for detailed information on the reagents.

Preparing Reagents

Dilute SPW Buffer with 100% ethanol as follows and store at room temperature.

Kit	100% Ethanol to be Added
D2500-00	20 mL
D2500-01	100 mL
D2500-02	100 mL per bottle

Guideline for Vacuum Manifold

The following is required for use with the Vacuum/Spin Protocol:

A) Vacuum Manifold

Compatible Vacuum Manifolds: Qiagen QIAvac24, Sigma Aldrich VM20, Promega Vacman®, or manifold with standard Luer connector


B) Vacuum Flask

C) Vacuum Tubing

D) Vacuum Source (review tables below for pressure settings)

Conversion from millibars:	Multiply by:
Millimeters of mercury (mm Hg)	0.75
Kilopascals (kPa)	0.1
Inches of mercury (inch Hg)	0.0295
Torrs (Torr)	0.75
Atmospheres (atmos)	0.000987
Pounds per Square Inch (psi)	0.0145

Illustrated Vacuum Setup:



Spin Protocol

Materials and Equipment to be Supplied by User:

- Heat block or water bath capable of 60°C
- Microcentrifuge capable of at least 13,000g
- Vortexer
- Nuclease-free 1.5 mL microcentrifuge tubes
- 100% ethanol
- Optional: 5M Sodium Acetate, pH 5.2
- Optional: Sterile deionized water

Before starting:

- Prepare SPW Buffer according to the “Preparing Reagents” section on Page 6
- Set heating block or water bath to 60°C

Note: The yellow color of the XP2 Binding Buffer signifies a pH of ≤ 7.5 .

1. Perform agarose gel/ethidium bromide electrophoresis to fractionate DNA fragments. Any type or grade of agarose may be used. However, it is strongly recommended that fresh TAE buffer or TBE buffer be used as running buffer. Do not reuse running buffer as its pH will increase and reduce yields.
2. When adequate separation of bands has occurred, carefully excise the DNA fragment of interest using a wide, clean, sharp scalpel. Minimize the size of the gel slice by removing extra agarose.
3. Determine the appropriate volume of the gel slice by weighing it in a clean 1.5 mL microcentrifuge tube. Assuming a density of 1 g/mL, the volume of gel is derived as follows: a gel slice of mass 0.3 g will have a volume of 0.3 mL.
4. Add 1 volume XP2 Binding Buffer.
5. Incubate at 60°C for 7 minutes or until the gel has completely melted. Vortex or shake the tube every 2-3 minutes.

E.Z.N.A.® Gel Extraction Kit

Important: Monitor the pH of the Gel/Binding Buffer mixture after the gel has completely dissolved. DNA yields will significantly decrease when the pH > 8.0. If the color of the mixture becomes orange or red, add 5 µL 5M sodium acetate (pH 5.2) to bring the pH down. After this adjustment, the color of the Gel/Binding Buffer mixture should be light yellow.

6. Insert a HiBind® DNA Mini Column in a 2 mL Collection Tube.
7. Add no more than 700 µL DNA/agarose solution from Step 5 to the HiBind® DNA Mini Column.
8. Centrifuge at 10,000*g* for 1 minute at room temperature.
9. Discard the filtrate and reuse collection tube.
10. Repeat Steps 7-9 until all of the sample has been transferred to the column.
11. Add 300 µL XP2 Binding Buffer.
12. Centrifuge at maximum speed ($\geq 13,000g$) for 1 minute at room temperature.
13. Discard the filtrate and reuse collection tube.
14. Add 700 µL SPW Buffer.

Note: SPW Buffer must be diluted with 100% ethanol prior to use. Please see Page 6 for instructions.
15. Centrifuge at maximum speed for 1 minute at room temperature.
16. Discard the filtrate and reuse collection tube.

E.Z.N.A.® Gel Extraction Kit

Optional: Repeat Steps 14-16 for a second SPW Buffer wash step. Perform the second wash step for any salt sensitive downstream applications.

17. Centrifuge the empty HiBind® DNA Mini Column for 2 minutes at maximum speed to dry the column matrix.

Note: It is important to dry the HiBind® DNA Mini Column matrix before elution. Residual ethanol may interfere with downstream applications.

18. Transfer the HiBind® DNA Mini Column to a clean 1.5 mL microcentrifuge tube.

19. Add 30-50 µL Elution Buffer or deionized water directly to the center of the column membrane.

Note: The efficiency of eluting DNA from the HiBind® DNA Mini Column is dependent on pH. If eluting DNA with deionized water, make sure that the pH is around 8.5.

20. Let sit at room temperature for 2 minutes.

21. Centrifuge at maximum speed for 1 minute.

Note: This represents approximately 70% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.

22. Store DNA at -20°C.

Vacuum Protocol

Materials and Equipment to be Supplied by User:

- Vacuum Manifold (See Page 7)
- Heat block or water bath capable of 60°C
- Microcentrifuge capable of at least 13,000g
- Vortexer
- Nuclease-free 1.5 mL microcentrifuge tubes
- 100% ethanol
- Optional: 5M Sodium Acetate, pH 5.2
- Optional: Sterile deionized water

Before starting:

- Prepare the Vacuum Manifold (See Page 7)
- Prepare SPW Buffer according to the “Preparing Reagents” section on Page 6
- Set heating block or water bath to 60°C

Note: The yellow color of the XP2 Binding Buffer signifies a pH of ≤ 7.5 .

1. Perform agarose gel/ethidium bromide electrophoresis to fractionate DNA fragments. Any type or grade of agarose may be used. However, it is strongly recommended that fresh TAE buffer or TBE buffer be used as running buffer. Do not reuse running buffer as its pH will increase and reduce yields.
2. When adequate separation of bands has occurred, carefully excise the DNA fragment of interest using a wide, clean, sharp scalpel. Minimize the size of the gel slice by removing extra agarose.
3. Determine the appropriate volume of the gel slice by weighing it in a clean 1.5 mL microcentrifuge tube. Assuming a density of 1 g/mL, the volume of gel is derived as follows: a gel slice of mass 0.3 g will have a volume of 0.3 mL.
4. Add 1 volume XP2 Binding Buffer.

E.Z.N.A.® Gel Extraction Kit

5. Incubate at 60°C for 7 minutes or until the gel has completely melted. Vortex or shake the tube every 2-3 minutes.
6. Prepare the vacuum manifold according to manufacturer's instructions.
7. Connect the HiBind® DNA Mini Column to the vacuum manifold.
8. Add no more than 700 µL DNA/agarose solution from Step 5 to the HiBind® DNA Mini Column.
9. Turn on the vacuum source to draw the sample through the column.
10. Turn off the vacuum.
11. Repeat Steps 8-10 until all of the sample has been transferred to the column.
12. Add 300 µL XP2 Binding Buffer.
13. Turn on the vacuum source to draw the sample through the column.
14. Turn off the vacuum.
15. Add 700 µL SPW Buffer.

Note: SPW Buffer must be diluted with 100% ethanol prior to use. Please see Page 6 for instructions.
16. Turn on the vacuum source to draw the sample through the column.
17. Turn off the vacuum.

E.Z.N.A.® Gel Extraction Kit

Optional: Repeat Steps 15-17 for a second SPW Buffer wash step. Perform the second wash step for any salt sensitive downstream applications.

18. Transfer the HiBind® DNA Mini Column to a clean 1.5 mL microcentrifuge tube.
19. Centrifuge the empty HiBind® DNA Mini Column for 2 minutes at maximum speed to dry the column matrix.

Note: It is important to dry the HiBind® DNA Mini Column matrix before elution. Residual ethanol may interfere with downstream applications.

20. Transfer the HiBind® DNA Mini Column to a clean 1.5 mL microcentrifuge tube.
21. Add 30-50 µL Elution Buffer or deionized water directly to the center of the column membrane.

Note: The efficiency of eluting DNA from the HiBind® DNA Mini Column is dependent on pH. If eluting DNA with deionized water, make sure that the pH is around 8.5.

22. Let sit at room temperature for 2 minutes.
23. Centrifuge at maximum speed for 1 minute.

Note: This represents approximately 70% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.

24. Store DNA at -20°C.

E.Z.N.A.[®] Gel Extraction Kit

Purification of DNA from Enzymatic Reactions

The following protocol is designed for DNA recovery from enzymatic reactions such as PCR and probe labeling reactions.

Materials and Equipment to be Supplied by User:

- Microcentrifuge capable of at least 13,000g
- Vortexer
- Nuclease-free 1.5 mL microcentrifuge tubes
- 100% ethanol

Before starting:

- Prepare SPW Buffer according to the “Preparing Reagents” section on Page 6

Note: The yellow color of the XP2 Binding Buffer signifies a pH of ≤ 7.5 .

1. Determine the volume of the enzymatic reaction.
2. Transfer the sample into a clean 1.5 mL microcentrifuge tube.
3. Add 1 volume XP2 Binding Buffer.
4. Vortex or invert the sample to mix thoroughly.
5. Briefly centrifuge the tube to collect any drops from the inside of the lid.
6. Insert a HiBind[®] DNA Mini Column into a 2 mL Collection Tube.
7. Transfer the sample to the HiBind[®] DNA Mini Column.
8. Centrifuge at 10,000g for 1 minute at room temperature.

E.Z.N.A.® Gel Extraction Kit

9. Discard the filtrate and reuse the collection tube.

10. Add 700 μ L SPW Buffer.

Note: SPW Buffer must be diluted with 100% ethanol prior to use. Please see Page 6 for instructions.

11. Centrifuge at maximum speed ($\geq 13,000g$) for 1 minute at room temperature.

12. Discard the filtrate and reuse collection tube.

13. Repeat Steps 10-12 for a second SPW Buffer wash step.

14. Centrifuge the empty HiBind® DNA Mini Column for 2 minutes at maximum speed to dry the column matrix.

Note: It is important to dry the HiBind® DNA Mini Column matrix before elution. Residual ethanol may interfere with downstream applications.

15. Transfer the HiBind® DNA Mini Column to a clean 1.5 mL microcentrifuge tube.

16. Add 30-50 μ L Elution Buffer or deionized water directly to the center of the column membrane.

Note: The efficiency of eluting DNA from the HiBind® DNA Mini Column is dependent on pH. If eluting DNA with deionized water, make sure that the pH is around 8.5.

17. Let sit at room temperature for 2 minutes.

18. Centrifuge at maximum speed for 1 minute.

Note: This represents approximately 70% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.

19. Store DNA at -20°C .

Troubleshooting Guide

Please use this guide to troubleshoot any problems that may arise. For further assistance, please contact the technical support staff, toll free, at (800-832-8896).

Possible Problems and Suggestions

Low DNA Yields	
Too little XP2 Binding Buffer added to gel	Volume of agarose gel slice determined incorrectly. Add enough XP2 Binding Buffer as instructed.
Agarose gel does not completely dissolve	Make sure water bath is set to 60°C and allow gel to completely melt. Add more XP2 Binding Buffer if necessary.
Inappropriate elution buffer	Use Elution Buffer or check the pH of the water used to elute DNA.
TAE/TBE running buffer is not fresh	With overuse, TAE buffer loses its buffering capacity and its pH increases. This raises the pH of the agarose/ DNA Binding Buffer solution which interferes with DNA binding to the HiBind® matrix. Adjust pH by adding 5 µL 5M sodium acetate, pH 5.2, to the gel slice. Use freshly prepared TAE buffer for gel purification in order to prevent the contamination of isolated DNA and improve yields.
Column clogged	
Agarose gel not completely dissolved in XP2 Binding Buffer	Make sure water bath is set to 55-60°C and allow gel to completely melt. For large agarose slices (>0.3 mL), it is recommended that the gel be diced into smaller fragments to aid melting.
No DNA eluted	
SPW Buffer not diluted with 100% ethanol	Prepare SPW Buffer as instructed on page 6 or as indicated on bottle.
Incorrect amount of Binding Buffer added	Measure the gel accurately and use 0.1 mL Binding Buffer (XP2) per 0.1 g gel.
Optical densities do not agree with DNA yield on agarose gel	
Trace contaminants eluted from column increase A_{260}	Make sure to wash column as instructed in Step 14 of the spin protocol, and Step 15 of the vacuum protocol. Alternatively, rely on agarose gel/ethidium bromide electrophoresis for quantization.
DNA sample floats out of well while loading agarose gel.	
Ethanol not completely removed from column	Centrifuge as instructed in Step 17 of the spin protocol and Step 19 of the vacuum protocol.

For more purification solutions, visit www.omegabiotek.com

AVAILABLE FORMATS



Spin Columns



96-Well
Silica Plates



Mag Beads

SAMPLE TYPES



Blood / Plasma



Plasmid



Cultured Cells



Plant & Soil



NGS Clean Up



Tissue



FFPE



Fecal Matter



BIO-TEK

innovations in nucleic acid isolation



Omega Bio-tek, Inc.
400 Pinnacle Way, Suite 450
Norcross, GA 30071



www.omegabiotek.com



770-931-8400



770-931-0230



info@omegabiotek.com



[omega-bio-tek](https://www.linkedin.com/company/omega-bio-tek)



[omegabiotek](https://twitter.com/omegabiotek)



[omegabiotek](https://www.facebook.com/omegabiotek)