



Mag-Bind® FFPE DNA 96 Kit

M6958-00	1 x 96 preps
M6958-01	4 x 96 preps

Manual Date: April 2022
Revision Number: v7.2

For Research Use Only

Mag-Bind® FFPE DNA 96 Kit

Table of Contents

Introduction	2
Kit Contents/Storage and Stability.....	3
Preparing Reagents.....	4
Mag-Bind® FFPE DNA 96 Protocol.....	5
Mag-Bind® FFPE DNA 96 Protocol with Xylene.....	9
Troubleshooting Guide.....	13
Notices & Disclaimers	14

Manual Date: April 2022

Revision Number: v7.2



Introduction

The Mag-Bind® FFPE DNA 96 kits provide a rapid and easy method for the isolation of total DNA from formalin-fixed, paraffin-embedded (FFPE) tissue sections. Due to fixation and embedding procedures, nucleic acids in FFPE samples are heavily fragmented and modified by formaldehyde. While the Mag-Bind® FFPE DNA 96 kits are optimized to minimize the effect of the formaldehyde modification, it is not recommended to use the DNA purified with these kits for downstream applications that require full length DNA.

The Mag-Bind® FFPE DNA 96 kits combine the high efficiency binding properties of Mag-Bind® technology with a specially designed buffer system to isolate total DNA sample from FFPE samples. There are two protocols included in this manual. The standard protocol uses a heating step to remove paraffin from the sample. The alternative protocol uses the traditional xylene extraction to remove paraffin.

After the paraffin removal steps, samples are first lysed in FTL2 Buffer with digestion of Proteinase K. The lysate is then heated to denature the proteinase and mixed with MB3 Buffer and magnetic particles to bind the nucleic acid on the surface of the Mag-Bind® particles. After two wash steps, purified DNA is eluted with Elution Buffer or nuclease-free water.

Important:

1. If automating this procedure on a liquid handler or a magnetic processor, please contact your Omega Bio-tek representative for instrument-specific instructions.
2. Kits include enough reagents for the specified number of preparations plus an additional 10% overage to ensure there is sufficient volume. Please be aware that the actual number of preparations may be lower due to pre-aliquoting of reagents, processing partial plates, and automation platform used etc. Additional reagents are available for purchase separately. Please visit the product page at www.omegabiotek.com or contact your Omega Bio-tek representative for more details and ordering information.

New in this Edition:

April 2022

- An important statement is included clarifying how the actual number of preparations is dependent on various factors and may be lower than the number of preparations specified with the kit.

January 2018

- RNase A is no longer supplied with this kit. An optional RNase A digestion step has been added.

Kit Contents

Mag-Bind® FFPE DNA 96 Kit	M6958-00	M6958-01
Preparations	1 x 96 preps	4 x 96 preps
Mag-Bind® Particles CH	3.3 mL	14 mL
FTL2 Buffer	30 mL	110 mL
MB4 Buffer	75 mL	250 mL
MPW Buffer	30 mL	125 mL
DNA Wash Buffer	50 mL	200 mL
LPA	1.1 mL	4.4 mL
Proteinase K Solution	2.2 mL	9 mL
Elution Buffer	30 mL	60 mL
User Manual	✓	✓

Storage and Stability

All of the Mag-Bind® FFPE DNA 96 Kit components are guaranteed for at least 12 months from the date of purchase when stored as follows. Mag-Bind® Particles CH should be stored at 2-8°C for long-term use. Proteinase K Solution can be stored at room temperature for up to 12 months. For long-term storage, store Proteinase K Solution at 2-8°C. All remaining components should be stored at room temperature. During shipment or storage in cool ambient conditions, precipitates may form. Dissolve such deposits by warming the solution at 37°C and gently shaking

Preparing Reagents

Dilute DNA Wash Buffer with 100% ethanol as follows and store at room temperature.

Kit	100% Ethanol to be Added
M6958-00	200 mL
M6958-01	800 mL

Dilute MPW Buffer with 100% isopropanol as follows and store at room temperature.

Kit	100% Isopropanol to be Added
M6958-00	30 mL
M6958-01	125 mL

Mag-Bind® FFPE DNA Protocol - 96-well plate without xylene

Important: If automating this procedure on a liquid handler or a magnetic processor, please contact your Omega Bio-tek representative for instrument-specific instructions.

Materials and Equipment to be Supplied by User:

- Microcentrifuge capable of $\geq 4,000g$
- Magnetic separation device for 96-well plates
- Water bath or heat block capable of 70°C
- Water bath or heat block capable of 80°C
- Vortexer
- 1.2 mL or 2.0 mL round-well plate
- Nuclease-free 96-well microplate
- 100% isopropanol
- 100% ethanol

Before Starting:

- Prepare Buffers according to Preparing Reagents section on Page 3.
- Set water baths or heat blocks to , 80°C, and 70°C.
- Vortex the Mag-Bind® Particles CH thoroughly before use.

1. Add 250 μ L FTL2 Buffer into a 1.2 mL or 2.0 mL round-well plate (not provided).
2. Cut 3-8 paraffin sample sections between 5-10 μ m.

Note: Do not use the first 2-3 sections from the sample block.

3. Immediately add 2-5 sections to the FTL2 Buffer.
4. Centrifuge at maximum speed ($\geq 4,000g$) for 5 minutes.
5. Incubate at 80°C for 15 minutes. Mix the sample a few times by gently vortexing the plate for 15 seconds. Make sure that the tissue sections stay submerged in the solution.
6. Let sit at room temperature for 5 minutes.

Mag-Bind® FFPE DNA Protocol - 96-well Plate without Xylene

7. Add 20 μ L Proteinase K Solution.
8. Incubate at 70°C for 3-5 hours with occasional mixing. If necessary, extend the incubation to overnight or until the tissue is completely lysed.

Optional: If RNA-free genomic DNA is required, add 10 μ L RNase A (20 mg/mL, not provided) and let sit for 5 minutes at room temperature.

9. Centrifuge at maximum speed ($>4,000g$) for 5 minutes. The paraffin will form a thin layer on top of the lysate solution.

10. Transfer 200 μ L cleared lysate into a new 1.2 mL or 2.0 mL round-well plate.

Tip: Use a 1 mL pipette tip or large orifice tip to penetrate the paraffin layer.

11. Add 500 μ L MB4 and 30 μ L Mag-Bind® Particles CH. Mix thoroughly by vortexing or pipetting up and down 10-20 times.
12. Let sit at room temperature for 5-10 minutes.
13. Place the plate on a magnetic separation device to magnetize the Mag-Bind® Particles CH. Let sit at room temperature until the Mag-Bind® Particles CH are completely cleared from solution.
14. Aspirate and discard the cleared supernatant. Do not disturb the Mag-Bind® Particles CH.
15. Remove the plate containing the Mag-Bind® Particles CH from the magnetic separation device.
16. Add 400 μ L MPW Buffer. Resuspend the Mag-Bind® Particles CH by vortexing or pipetting up and down 20 times.

Note: MPW Buffer must be diluted with 100% isopropanol prior to use. Please see Page 4 for instructions.

Mag-Bind® FFPE DNA Protocol - 96-well Plate without Xylene

17. Place the plate on a magnetic separation device to magnetize the Mag-Bind® Particles CH. Let sit at room temperature until the Mag-Bind® Particles CH are completely cleared from solution.
18. Aspirate and discard the cleared supernatant. Remove any liquid drops from each well. Do not disturb the Mag-Bind® Particles CH.
19. Remove the plate containing the Mag-Bind® Particles CH from the magnetic separation device.
20. Add 400 µL DNA Wash Buffer. Resuspend the Mag-Bind® Particles CH by vortexing or pipetting up and down 20 times.

Note: DNA Wash Buffer must be diluted with ethanol prior to use. Please see Page 4 for instructions.

21. Place the plate on a magnetic separation device to magnetize the Mag-Bind® Particles CH. Let sit at room temperature until the Mag-Bind® Particles CH are completely cleared from solution.
 22. Aspirate and discard the cleared supernatant. Do not disturb the Mag-Bind® Particles CH.
 23. Repeat Steps 20-22 for a second DNA Wash step.
 24. Select one of the following ethanol removal steps:
 - A. Leave the plate on the magnetic separation device. Add 400 µL nuclease-free water (not provided), leave on magnet for 20-30 seconds, and then aspirate. Do not leave nuclease-free water on Mag-Bind® Particles CH for more than 60 seconds. Continue to Step 25.
- OR**
- B. Leave the plate on the magnetic separation device. Wait 1 minute. Remove residual liquid with a pipettor. Dry the Mag-Bind® Particles CH for an additional 10 minutes. Continue to Step 25.

Mag-Bind® FFPE DNA Protocol - 96-well Plate without Xylene

25. Remove the plate containing the Mag-Bind® Particles CH from the magnetic separation device.
26. Add 30-50 μ L Elution Buffer. Resuspend the Mag-Bind® Particles CH by vortexing or pipetting up and down 30 times.
27. Incubate at room temperature for 5 minutes.
28. Place the plate on a magnetic separation device to magnetize the Mag-Bind® Particles CH. Let sit at room temperature until the Mag-Bind® Particles CH are completely cleared from solution.
29. Transfer the cleared supernatant containing purified DNA to a nuclease-free 96-well microplate (not provided). Store the DNA at -20°C.

Mag-Bind® FFPE DNA Protocol - 96-well Plate with Xylene

Mag-Bind® FFPE DNA 96 - 96-well plate with xylene

Important: If automating this procedure on a liquid handler or a magnetic processor, please contact your Omega Bio-tek representative for instrument-specific instructions.

Note: The following protocol uses xylene to remove paraffin from the FFPE sample.
Use fume hood and take proper protection during xylene extraction.

Materials and Equipment to be Supplied by User:

- Centrifuge with swing-bucket rotor capable of 4,000g
- Rotor adaptor for 96-well deep-well plates
- Magnetic separation device for 96-well deep-well plates
- Water bath or heat block capable of 55°C
- Water bath or heat block capable of 80°C
- Water bath or heat block capable of 90°C
- Vortexer
- 1.2 mL or 2.0 mL round-well plates
- Nuclease-free 96-well microplates
- 100% isopropanol
- 100% ethanol
- Xylene

Before Starting:

- Prepare Buffers according to Preparing Reagents section on Page 3.
- Set water baths or heat blocks to 90°C, 80°C, and 55°C.
- Vortex the Mag-Bind® Particles CH thoroughly before use.

1. Add 1 mL xylene (not provided) into each well of a 1.2 mL or 2.0 mL round-well plate (not provided).
2. Cut 3-8 paraffin sample sections between 5-10 µm.

Note: Do not use the first 2-3 sections from the sample block.

3. Immediately add 2-5 sections to the xylene. Vortex for 20 seconds to mix thoroughly.

Mag-Bind® FFPE DNA Protocol - 96-well Plate with Xylene

4. Centrifuge at 4,000g for 5-10 minutes to pellet the tissue.

Note: If the tissue does not form a tight pellet, centrifuge for an additional 5 minutes.

5. Aspirate and discard the xylene. Do not disturb the tissue pellet.
6. Add 1 mL 100% ethanol (not provided). Vortex for 20 seconds to mix thoroughly.
7. Centrifuge at 4,000g for 5 minutes to pellet the tissue. The pellet should appear opaque.
8. Aspirate and discard the ethanol. Do not disturb the tissue pellet. Remove any liquid drops with a pipette.
9. Repeat Steps 6-8 for a second ethanol wash step.
10. Let sit at room temperature for 10-20 minutes.

Note: It is critical to completely dry the sample before the next Proteinase K digestion step. Residual ethanol will affect the efficiency of the Proteinase K digestion. If a vacuum oven is available, place the plate in the vacuum oven preset at 45°C for 10-20 minutes.
11. Add 250 µL FTL2 Buffer and 20 µL Proteinase K Solution. Resuspend the pellet by vortexing or pipetting up and down 20 times.
12. Incubate at 55°C for 3-5 hours with occasional mixing. If necessary, extend the incubation to overnight or until the tissue is completely lysed.
13. Incubate at 90°C for 45-60 minutes.

Optional: If RNA-free genomic DNA is required, add 10 µL RNase A (20 mg/mL, not provided) and let sit for 5 minutes at room temperature.
14. Centrifuge at 4,000g for 5 minutes.

Mag-Bind® FFPE DNA Protocol - 96-well Plate with Xylene

15. Transfer 200 μ L cleared supernatant into a new 1.2 mL or 2.0 mL round-well plate.
16. Add 500 μ L MB4 Buffer and 30 μ L Mag-Bind® Particles CH. Mix thoroughly by vortexing or pipetting up and down 10-20 times.

Note: If DNA content from sample is expected to be low, add 10 μ L LPA.

17. Let sit at room temperature for 5-10 minutes.
18. Place the plate on the magnetic separation device to magnetize the Mag-Bind® Particles CH. Let sit at room temperature until the Mag-Bind® Particles CH are completely cleared from solution.
19. Aspirate and discard the supernatant. Do not disturb the Mag-Bind® Particles CH.
20. Remove the plate from the magnetic separation device.
21. Add 400 μ L MPW Buffer. Resuspend the Mag-Bind® Particles CH by vortexing or pipetting up and down 20 times.

Note: MPW Buffer must be diluted with 100% isopropanol prior to use. Please see Page 4 for instructions.

22. Place the plate on the magnetic separation device to magnetize the Mag-Bind® Particles CH. Let sit at room temperature until the Mag-Bind® Particles CH are completely cleared from solution.
23. Aspirate and discard the supernatant. Do not disturb the Mag-Bind® Particles CH.
24. Remove the plate from the magnetic separation device.
25. Add 400 μ L DNA Wash Buffer. Resuspend the Mag-Bind® Particles CH by vortexing or pipetting up and down 20 times.

Note: DNA Wash Buffer must be diluted with ethanol prior to use. Please see Page 4 for instructions.

Mag-Bind® FFPE DNA Protocol - 96-well Plate with Xylene

26. Place the plate on the magnetic separation device to magnetize the Mag-Bind® Particles CH. Let sit at room temperature until the Mag-Bind® Particles CH are completely cleared from solution.
27. Aspirate and discard the cleared supernatant. Do not disturb the Mag-Bind® Particles CH.
28. Repeat Steps 25-27 for a second DNA Wash step.
29. Select one of the following ethanol removal steps:
 - A. Leave the plate on the magnetic separation device. Add 400 µL nuclease-free water (not provided), leave on magnet for 20-30 seconds, and then aspirate. Do not leave nuclease-free water on Mag-Bind® Particles CH for more than 60 seconds. Continue to Step 30.

OR

- B. Leave the plate on the magnetic separation device. Wait 1 minute. Remove residual liquid with a pipettor. Dry the Mag-Bind® Particles CH for an additional 10 minutes. Continue to Step 30.
30. Remove the plate containing the Mag-Bind® Particles CH from the magnetic separation device.
31. Add 30-50 µL Elution Buffer. Resuspend the Mag-Bind® Particles CH by vortexing or pipetting up and down 30 times.
32. Incubate at room temperature for 5 minutes.
33. Place the plate on a magnetic separation device to magnetize the Mag-Bind® Particles CH. Let sit at room temperature until the Mag-Bind® Particles CH are completely cleared from solution.
34. Transfer the cleared supernatant containing purified DNA to a nuclease-free 96-well microplate (not provided). Store the DNA at -20°C.

Troubleshooting Guide

Please use this guide to troubleshoot any problems that may arise. For further assistance, please contact the technical support staff, toll free, at **1-800-832-8896**.

Problem	Cause	Solution
Low DNA yields	Incomplete resuspension of magnetic particles	Resuspend the magnetic particles by vortexing before use.
	DNA degraded during sample storage	Make sure the sample is properly stored and make sure the samples are processed immediately after collection or removal from storage.
	MPW Buffer and DNA Wash Buffer were not prepared correctly	Prepare MPW Buffer and DNA Wash Buffer according to the instructions on Page 4.
	Loss of magnetic beads during operation	Increase the bead collection time.
Problem	Cause	Solution
Problem with downstream application	DNA is over fixated during tissue formalin fixation	Extend incubation time at 90°C to 90 minutes.
Carryover of the magnetic beads in the elution	Carryover the magnetic beads in the eluted DNA will not effect downstream applications	To remove the carryover magnetic particles from the eluted DNA, simply magnetize the magnetic particles and carefully transfer the DNA eluate to a new plate.

Notices & Disclaimers

For European Union Use.

MB4 Buffer contains Triton X-100, 2-[4-(2,4,4-trimethylpentan-2-yl)phenoxy]ethanol (CAS 9002-93-1), a substance included in the European Authorisation list (Annex XIV) of REACH Regulation (EC) No 1907/2006. Substances and mixtures used for the purpose of Scientific Research and Development (SR&D) are exempt from authorization requirements if used below 1 tonne per year in volume.

Scientific Research and Development includes experimental research or analytical activities at a laboratory scale such as synthesis and testing of applications of chemicals, release tests, etc. as well as the use of the substance in monitoring and routine quality control or in vitro diagnostics.

Notes:

Notes:

For more purification solutions, visit www.omegabiotek.com

AVAILABLE FORMATS



Spin Columns



96-Well
Silica Plates



Mag Beads

SAMPLE TYPES



Blood / Plasma



Plasmid



Cultured Cells



Plant & Soil



NGS Clean Up



Tissue



FFPE



Fecal Matter



BIO-TEK

innovations in nucleic acid isolation



Omega Bio-tek, Inc.
400 Pinnacle Way, Suite 450
Norcross, GA 30071



www.omegabiotek.com



770-931-8400



770-931-0230



info@omegabiotek.com



[omega-bio-tek](https://www.linkedin.com/company/omega-bio-tek)



[omegabiotek](https://twitter.com/omegabiotek)



[omegabiotek](https://www.facebook.com/omegabiotek)