



E.Z.N.A.[®] Blood DNA Mini Kit

D3392-02	200 preps
D3392-03	600 preps

Manual Date: April 2024
Revision Number: v8.2

For Research Use Only

E.Z.N.A.® Blood DNA Mini Kit

Table of Contents

Introduction	2
Illustrated Protocol.....	3
Kit Contents/Storage and Stability.....	4
Warnings and Safety Information.....	5
Preparing Reagents.....	6
Blood and Body Fluid Protocol.....	7
Buccal Swab Protocol.....	10
Dried Blood Protocol.....	13
Buffy Coat Protocol.....	16
Troubleshooting Guide.....	19

Manual Date: April 2024

Revision Number: v8.2



Introduction

The E.Z.N.A.® family of products is an innovative system that radically simplifies the extraction and purification of nucleic acids from a variety of sources. The key to this system is the new HiBind® matrix that specifically, but reversibly, binds DNA or RNA under optimal conditions allowing proteins and other contaminants to be removed. Nucleic acids are easily eluted with deionized water or a low salt buffer.

The E.Z.N.A.® Blood DNA Mini Kit provides an easy and rapid method for the isolation of genomic DNA for consistent PCR and Southern analysis. Up to 250 µL fresh, frozen, or anticoagulated whole blood can be readily processed at one time. The E.Z.N.A.® Blood DNA Mini Kit can also be used for the preparation of genomic DNA from buffy coat, serum, plasma, saliva, buccal swabs, and other body fluids. The E.Z.N.A.® Blood DNA Kit allows for single or multiple simultaneous processing of multiple samples. There is no need for phenol/chloroform extractions, and time-consuming steps are eliminated (e.g. precipitation using isopropanol or ethanol). Purified DNA obtained with the E.Z.N.A.® Blood DNA Kit is ready for applications such as PCR, restriction digestion, and Southern blotting.

Benefits of the E.Z.N.A.® Blood DNA Mini Kit

- Optimized buffers that guarantee pure DNA
- No organic extractions
- Purified DNA can be directly used in most downstream applications

Binding Capacity

Each HiBind® DNA Mini Column can bind approximately 100 µg DNA. Using greater than 250 µL whole blood or buffy coat is not recommended.

New in this Edition:

April 2024

- Addition of Warnings and Safety Information.

June 2023:

- General revision.

March 2022:

- A new 600 prep kit (D3392-03) has been added and is now available for purchase.


Illustrated Protocol



Lyse



Adjust Binding
Conditions



Bind



Wash 3X



Dry



Elute

Kit Contents

Product	D3392-02	D3392-03
Purifications	200 preps	600 preps
HiBind® DNA Mini Columns	200	600
2 mL Collection Tubes	400	1200
BL Buffer	60 mL	3 x 110 mL
HBC Buffer	80 mL	250 mL
DNA Wash Buffer	3 x 25 mL	200 mL
Elution Buffer	250 mL	300 mL
Proteinase K Solution	6 mL	18 mL
User Manual	✓	✓

Storage and Stability

All of the E.Z.N.A.® Blood DNA Mini Kit components are guaranteed for 12 months from the date of purchase when stored as follows. Proteinase K Solution can be stored at room temperature for 12 months. For long-term storage (>12 months), store Proteinase K Solution at 2-8°C. Store all other components at room temperature and away from bright light. During shipment or storage under cool ambient conditions, a precipitate may form in some buffers. Dissolve such deposits by warming the solution to 37°C and gently shaking.

Warnings and Safety Information

Warnings

This kit is for research use only.

Please read all instructions carefully before using the kit.

Decontaminate and dispose of all potentially infectious materials in accordance with applicable local, state, and national regulations. Please refer to safety data sheets (SDSs) for information on disposal of different components included in this kit.

Safety Information

All chemicals and biological materials are potentially hazardous. Biological samples such as plasma, serum, tissues, body fluids, blood etc. are potentially infectious and must be treated as biohazardous materials. Conduct all work in properly equipped facilities following universal precautions and using appropriate personal safety equipment such as disposable gloves, lab coats, safety glasses etc. as required by policies and procedures outlined by your facility. Please refer to safety data sheets (SDSs) for information on safe handling, transport and disposal of different components included in this kit. SDSs are made available in PDF format on the product page at www.omegabiotek.com. Discard all waste in accordance with the local safety regulations.

Some of the buffers included in the product contain guanidine-based chaotropic agents, which can form highly reactive compounds when combined with bleach. DO NOT add bleach or acidic solutions to guanidine-containing waste. Please access the SDSs online for detailed information on the reagents.

Preparing Reagents

1. Dilute HBC Buffer with 100% isopropanol as follows and store at room temperature.

Kit	100% Isopropanol to be Added
D3392-02	32 mL
D3392-03	100 mL

2. Dilute DNA Wash Buffer with 100% ethanol as follows and store at room temperature.

Kit	100% Ethanol to be Added
D3392-02	100 mL per bottle
D3392-03	800 mL

E.Z.N.A.® Blood DNA Mini Kit

Blood and Body Fluids Protocol

The procedure below has been optimized for the use with fresh or frozen blood samples up to 250 µL in volume. Anti-coagulated blood, saliva, serum, buffy coat, or other body fluids can also be used. In addition, $\leq 10^7$ of leukocytes or cultured cells may be used with this procedure. For DNA extraction from tissue and mouse tail we suggest that you use the E.Z.N.A.® Tissue DNA Kit (Product No. D3396). To isolate viral RNA from serum or other non-cellular body fluids the E.Z.N.A.® Viral RNA Kit (Product No. R6874) is recommended.

Materials and Equipment to be Supplied by User:

- Tabletop microcentrifuge capable of at least 13,000g
- Water bath, incubator, or heat block capable of 65°C
- Vortexer
- Nuclease-free 2 mL microcentrifuge tubes
- 100% ethanol
- 100% isopropanol
- Optional: 10mM Tris-HCl or PBS
- Optional: RNase stock solution (50 mg/mL; when RNA-free genomic DNA is required)

Before Starting:

- Prepare HBC Buffer and DNA Wash Buffer according to the directions in the "Preparing Reagents" section on Page 6
 - Set water bath, incubator, or heat block to 65°C
 - Heat the Elution Buffer to 65°C
-
1. Transfer the sample into a nuclease-free 2 mL microcentrifuge tube (not provided) and bring the volume up to 250 µL with 10mM Tris-HCl (not provided), PBS (not provided), or Elution Buffer.
 2. Add 25 µL Proteinase K Solution and 250 µL BL Buffer. Vortex at maximum speed for 15 seconds.

Optional: If RNA-free genomic DNA is required, add 5 µL RNase A (50 mg/mL).

3. Incubate for 10 minutes at 65°C. Vortex briefly once during incubation.
4. Add 260 µL 100% ethanol (not provided). Vortex at maximum speed for 20 seconds.

E.Z.N.A.[®] Blood DNA Mini Kit

5. Centrifuge briefly to collect any drops from the inside of the lid.
6. Insert a HiBind[®] DNA Mini Column into a 2 mL collection tube.
7. Transfer the entire sample to the column.
8. Centrifuge at $\geq 10,000g$ for 1 minute.
9. Discard the filtrate and the collection tube.
10. Insert the HiBind[®] DNA Mini Column into a new 2 mL collection tube.
11. Add 500 μ L HBC Buffer.

Note: HBC Buffer must be diluted with 100% isopropanol before use. Please see Page 6 for instructions.

12. Centrifuge at $\geq 10,000g$ for 1 minute.
13. Discard the filtrate and reuse the collection tube.
14. Add 700 μ L DNA Wash Buffer.

Note: DNA Wash Buffer must be diluted with 100% ethanol before use. Please see Page 6 for instructions.

15. Centrifuge at 10,000g for 1 minute.
16. Discard the filtrate and reuse the collection tube.

E.Z.N.A.® Blood DNA Mini Kit

17. Repeat Steps 14-16 for a second DNA Wash Buffer step.
18. Centrifuge the empty HiBind® DNA Mini Column at maximum speed ($\geq 10,000g$) for 2 minutes to dry the column matrix.

Note: It is important to dry the column membrane before elution. Residual ethanol may interfere with downstream applications.

19. Transfer the HiBind® DNA Mini Column into a nuclease-free 2 mL microcentrifuge tube.
20. Add 100-200 μ L Elution Buffer heated to 65°C.
21. Let sit for 5 minutes at room temperature.

Note: Incubating the HiBind® DNA Mini Column at 65°C rather than room temperature will give a modest increase in DNA yield per elution.

22. Centrifuge at $\geq 13,000g$ for 1 minute.
23. Repeat Steps 20-22 for a second elution step.

Note: Any combination of the following steps can be used to help increase DNA yield.

- After adding the Elution Buffer, incubate the column for 5 minutes.
- Increase the elution volume.
- Repeat the elution step with fresh Elution Buffer (this may increase the yield, but decrease the concentration).
- Repeat the elution step using the eluate from the first elution (this may increase yield while maintaining elution volume).

24. Store eluted DNA at -20°C.

E.Z.N.A.[®] Blood DNA Mini Kit

Buccal Swabs Protocol

This protocol requires an increased volume of BL Buffer. Fewer preparations can be performed. Additional BL Buffer can be purchased separately.

Materials and Equipment to be Supplied by User:

- Tabletop microcentrifuge capable of at least 13,000g
- Water bath, incubator, or heat block capable of 65°C
- Vortexer
- Nuclease-free 2 mL microcentrifuge tubes
- 100% ethanol
- 100% isopropanol
- PBS
- Optional: RNase stock solution (50 mg/mL; when RNA-free genomic DNA is required)

Before Starting:

- Prepare HBC Buffer and DNA Wash Buffer according to the directions in the “Preparing Reagents” section on Page 6
- Set water bath, incubator, or heat block to 65°C
- Heat the Elution Buffer to 65°C

1. Place the buccal swab in a nuclease-free 2 mL microcentrifuge tube (not provided).
2. Add 500 μ L PBS (not provided).

Optional: If RNA-free genomic DNA is required, add 5 μ L RNase A (50 mg/mL).

3. Add 25 μ L Proteinase K Solution and 500 μ L BL Buffer. Vortex at maximum speed for 30 seconds.
4. Incubate for 10 minutes at 65°C.
5. Discard the buccal swab.
6. Add 500 μ L 100% ethanol (not provided). Vortex at maximum speed for 20 seconds.

E.Z.N.A.® Blood DNA Mini Kit

7. Centrifuge briefly to collect any drops from the inside of the lid.
8. Insert a HiBind® DNA Mini Column into a 2 mL collection tube.
9. Transfer 750 µL sample to the column.
10. Centrifuge at 10,000g for 1 minute.
11. Discard the filtrate and reuse the collection tube.
12. Repeat Steps 9-11 until all the sample has been transferred to the column.
13. Place the HiBind DNA Mini Column into a new 2 mL collection tube.

14. Add 500 µL HBC Buffer.

Note: HBC Buffer must be diluted with 100% isopropanol before use. Please see Page 6 for instructions.

15. Centrifuge at $\geq 10,000g$ for 1 minute.
16. Discard the filtrate and reuse collection tube.
17. Add 700 µL DNA Wash Buffer.

Note: DNA Wash Buffer must be diluted with 100% ethanol before use. Please see Page 6 for instructions.

18. Centrifuge at 10,000g for 1 minute.
19. Discard the filtrate and reuse the collection tube.

E.Z.N.A.® Blood DNA Mini Kit

20. Repeat Steps 17-19 for a second DNA Wash Buffer step.
21. Centrifuge the empty HiBind® DNA Mini Column at maximum speed ($\geq 10,000g$) for 2 minutes to dry the column matrix.

Note: It is important to dry the column membrane before elution. Residual ethanol may interfere with downstream applications.

22. Transfer the HiBind® DNA Mini Column into a nuclease-free 2 mL microcentrifuge tube.
23. Add 100-200 μ L Elution Buffer heated to 65°C.

24. Let sit for 5 minutes at room temperature.

Note: Incubating the HiBind® DNA Mini Column at 65°C rather than room temperature will give a modest increase in DNA yield per elution.

25. Centrifuge at $\geq 13,000g$ for 1 minute.
26. Repeat Steps 23-25 for a second elution step.

Note: Any combination of the following steps can be used to help increase DNA yield.

- After adding the Elution Buffer, incubate the column for 5 minutes.
- Increase the elution volume.
- Repeat the elution step with fresh Elution Buffer (this may increase the yield, but decrease the concentration).
- Repeat the elution step using the eluate from the first elution (this may increase yield while maintaining elution volume).

27. Store eluted DNA at -20°C.

Dried Blood Protocol

Materials and Equipment to be Supplied by User:

- Tabletop microcentrifuge capable of at least 13,000g
- Water bath, incubator, or heat block capable of 65°C
- Vortexer
- Nuclease-free 2 mL microcentrifuge tubes
- 100% ethanol
- 100% isopropanol
- PBS
- Optional: RNase stock solution (50 mg/mL; when RNA-free genomic DNA is required)

Before Starting:

- Prepare HBC Buffer and DNA Wash Buffer according to the directions in the “Preparing Reagents” section on Page 6
 - Set water bath, incubator, or heat block to 65°C
 - Heat the Elution Buffer to 65°C
-
1. Cut or punch-out the blood spot from the filter paper (up to 200 µL blood can be used per spot). Tear or cut the filter paper into small pieces and place them into a nuclease-free 2 mL microcentrifuge tube (not provided).
 2. Add 250 µL PBS (not provided).
 3. Incubate for 1 hour at 65°C. Vortex briefly every 20 minutes.
 4. Add 25 µL Proteinase K Solution. Vortex at maximum speed for 15 seconds.
 5. Incubate for 30 minutes at 65°C. Vortex briefly several times during incubation.
 6. Centrifuge at $\geq 13,000g$ for 5 minutes.
 7. Transfer the supernatant to a nuclease-free 2 mL microcentrifuge tube.

E.Z.N.A.® Blood DNA Mini Kit

8. Add 1 volume BL Buffer and 1 volume 100% ethanol (not provided). Vortex to mix thoroughly.
9. Centrifuge briefly to collect any drops from the inside of the lid.
10. Insert a HiBind® DNA Mini Column into a 2 mL collection tube.
11. Transfer the entire sample to the column.
12. Centrifuge at $\geq 10,000g$ for 1 minute.
13. Discard the filtrate and the collection tube.
14. Insert the HiBind® DNA Mini Column into a new 2 mL collection tube.
15. Add 500 μ L HBC Buffer.

Note: HBC Buffer must be diluted with 100% isopropanol before use. Please see Page 6 for instructions.

16. Centrifuge at $\geq 10,000g$ for 1 minute.
17. Discard the filtrate and reuse the collection tube.
18. Add 700 μ L DNA Wash Buffer.

Note: DNA Wash Buffer must be diluted with 100% ethanol before use. Please see Page 6 for instructions.

19. Centrifuge at $10,000g$ for 1 minute.
20. Discard the filtrate and reuse the collection tube.

E.Z.N.A.® Blood DNA Mini Kit

21. Repeat Steps 18-20 for a second DNA Wash Buffer step.
22. Centrifuge the empty HiBind® DNA Mini Column at maximum speed ($\geq 10,000g$) for 2 minutes to dry the column matrix.

Note: It is important to dry the column membrane before elution. Residual ethanol may interfere with downstream applications.

23. Transfer the HiBind® DNA Mini Column into a nuclease-free 2 mL microcentrifuge tube.
24. Add 100-200 μ L Elution Buffer heated to 65°C.

25. Let sit for 5 minutes at room temperature.

Note: Incubating the HiBind® DNA Mini Column at 65°C rather than room temperature will give a modest increase in DNA yield per elution.

26. Centrifuge at $\geq 13,000g$ for 1 minute.
27. Repeat Steps 24-26 for a second elution step.

Note: Any combination of the following steps can be used to help increase DNA yield.

- After adding the Elution Buffer, incubate the column for 5 minutes.
- Increase the elution volume.
- Repeat the elution step with fresh Elution Buffer (this may increase the yield, but decrease the concentration).
- Repeat the elution step using the eluate from the first elution (this may increase yield while maintaining elution volume).

28. Store eluted DNA at -20°C.

Note: Blood spots from finger pricks usually contain no more than 50 μ L blood, and yield approximately 500 ng to 1 μ g DNA. This is usually sufficient for PCR analysis. To obtain higher DNA concentrations, elute with 50 μ L preheated elution buffer (volumes lower than 50 μ L greatly reduce yields). Alternatively, the first eluate can be used to perform a second elution.

E.Z.N.A.[®] Blood DNA Mini Kit

Buffy Coat Protocol

The buffy coat fraction of whole blood is enriched with leukocytes and usually gives at least 5-fold more DNA than the same volume of blood. To prepare the buffy coat from fresh whole blood, simply centrifuge the sample at 3,000-4,000g for 10 minutes at room temperature. Three layers should form: a plasma upper layer, a buffy coat middle layer, and an erythrocyte bottom layer. Carefully aspirate the plasma, making sure not to disturb the layer of concentrated leukocytes. The buffy coat can be drawn off with a pipette and used directly in the E.Z.N.A.[®] Blood DNA Mini Kit or frozen at -70°C.

This protocol requires an increased volume of BL Buffer. Fewer preparations can be performed. Additional BL Buffer can be purchased separately.

Materials and Equipment to be Supplied by User:

- Tabletop microcentrifuge capable of at least 13,000g
- Water bath, incubator, or heat block capable of 65°C
- Vortexer
- Nuclease-free 2 mL microcentrifuge tubes
- 100% ethanol
- 100% isopropanol
- Optional: 10mM Tris-HCl or PBS
- Optional: RNase stock solution (50 mg/mL; when RNA-free genomic DNA is required)

Before Starting:

- Prepare HBC Buffer and DNA Wash Buffer according to the directions in the “Preparing Reagents” section on Page 6
 - Set water bath, incubator, or heat block to 65°C
 - Heat the Elution Buffer to 65°C
1. Transfer the sample into a nuclease-free 2 mL microcentrifuge tube (not provided) and bring the volume up to 500 µL with 10mM Tris-HCl (not provided), PBS (not provided), or Elution Buffer.
 2. Add 25 µL Proteinase K Solution and 500 µL BL Buffer. Vortex at maximum speed for 15 seconds.

Optional: If RNA-free genomic DNA is required, add 2 µL RNase A (50 mg/mL).

3. Incubate for 10 minutes at 65°C. Vortex briefly once during incubation.

E.Z.N.A.® Blood DNA Mini Kit

4. Add 500 μ L 100% ethanol (not provided). Vortex at maximum speed for 20 seconds.
5. Centrifuge briefly to collect any drops from the inside of the lid.
6. Insert a HiBind® DNA Mini Column into a 2 mL collection tube.
7. Transfer the 750 μ L sample to the column.
8. Centrifuge at $\geq 10,000g$ for 1 minute.
9. Discard the filtrate and the collection tube.
10. Repeat Steps 7-9 until all the sample has been transferred to the column.
11. Insert the HiBind® DNA Mini Column into a new 2 mL collection tube.
12. Add 500 μ L HBC Buffer.

Note: HBC Buffer must be diluted with 100% isopropanol before use. Please see Page 6 for instructions.
13. Centrifuge at $\geq 10,000g$ for 1 minute.
14. Discard the filtrate and reuse the collection tube.
15. Add 700 μ L DNA Wash Buffer.

Note: DNA Wash Buffer must be diluted with 100% ethanol before use. Please see Page 6 for instructions.
16. Centrifuge at 10,000g for 1 minute.

E.Z.N.A.® Blood DNA Mini Kit

17. Discard the filtrate and reuse the collection tube.
18. Repeat Steps 15-17 for a second DNA Wash Buffer step.
19. Centrifuge the empty HiBind® DNA Mini Column at maximum speed ($\geq 10,000g$) for 2 minutes to dry the column matrix.

Note: It is important to dry the column membrane before elution. Residual ethanol may interfere with downstream applications.

20. Transfer the HiBind® DNA Mini Column into a nuclease-free 2 mL microcentrifuge tube.
21. Add 100-200 μ L Elution Buffer heated to 65°C.
22. Let sit for 5 minutes at room temperature.

Note: Incubating the HiBind® DNA Mini Column at 65°C rather than room temperature will give a modest increase in DNA yield per elution.

23. Centrifuge at $\geq 13,000g$ for 1 minute.
24. Repeat Steps 21-23 for a second elution step.

Note: Any combination of the following steps can be used to help increase DNA yield.

- After adding the Elution Buffer, incubate the column for 5 minutes.
- Increase the elution volume.
- Repeat the elution step with fresh Elution Buffer (this may increase the yield, but decrease the concentration).
- Repeat the elution step using the eluate from the first elution (this may increase yield while maintaining elution volume).

25. Store eluted DNA at -20°C.

Troubleshooting Guide

Please use this guide to troubleshoot any problems that may arise. For further assistance, please contact the technical support staff, toll free, at **1-800-832-8896**.

Problem	Cause	Solution
Clogged Column	Incomplete Lysis	Extend incubation time with BL Buffer and Proteinase K Solution.
	Too much Sample	Divide sample into multiple tubes and adjust the volume to 250 μ L with BL Buffer.
	Sample is Viscous	Divide sample into multiple tubes and adjust the volume to 250 μ L with BL Buffer.
Problem	Cause	Solution
Low DNA Yield	Poor Elution	Repeat elution with increased elution volume. Incubate columns at 65°C for 5 minutes with Elution Buffer.
	Improper Washing	DNA Wash Buffer must be diluted with 100% ethanol before use. If refrigerated, DNA Wash Buffer must be brought to room temperature.
		HBC Buffer must be diluted with isopropanol before use. If refrigerated, HBC Buffer must be brought to room temperature.
	Sample has low DNA Content	Increase starting material and volume of all reagents (Proteinase K Solution, BL Buffer, ethanol) proportionally. Load aliquots of lysate through column successively.
	Prime Columns	Add 100 μ L 3M NaOH to the column prior to loading the sample. Let sit for 4 minutes. Centrifuge at 10,000g for 60 seconds. Add 100 μ L water to the column and centrifuge at 10,000g for 60 seconds. Discard the filtrate.

Notes:

For more purification solutions, visit www.omegabiotek.com

AVAILABLE FORMATS



Spin Columns



96-Well
Silica Plates



Mag Beads

SAMPLE TYPES



Blood / Plasma



Plasmid



Cultured Cells



Plant & Soil



NGS Clean Up



Tissue



FFPE




Fecal Matter



BIO-TEK

innovations in nucleic acid isolation

 Omega Bio-tek, Inc.
400 Pinnacle Way, Suite 450
Norcross, GA 30071

 www.omegabiotek.com

 770-931-8400

 770-931-0230

 info@omegabiotek.com

 [omega-bio-tek](https://www.linkedin.com/company/omega-bio-tek)

 [omegabiotek](https://twitter.com/omegabiotek)

 [omegabiotek](https://www.facebook.com/omegabiotek)