

Revision Date: February 2024
Revision Number: v2.0

KWB Buffer

Product Number	KWB-450
KWB Buffer	450 mL
Instruction Manual	✓

Storage and Stability

KWB Buffer is stable for at least 12 months from date of purchase when stored at room temperature and away from bright light. Do not refrigerate. During shipment or storage in cool ambient conditions, precipitates may form. Dissolve such deposits by warming the solution at 37°C and gently shaking. Precipitates may also form when there is excess evaporation, so store container with the lid tightly closed.